Ortho Slice™ 3D Knee Planning

Security Operations Manual

REF: 6007-670-000


Table of Contents

01	Purpose	01
02	Definitions	01
	Product Description	04
3.1	Device and Manufacturer Identification	04
3.2	Device Intended Use	05
3.3	Vulnerability Intake and Monitoring	05
3.4	System Characterization and System Assets	05
3.5	System Security Context and Intended Environment	06
3.6	Setup of the SaMD (Software as a Medical Device) Ortho Slice 3D Knee planning software	07
04	Management of PII and PHI	07
4.1	Handling of Patient Requests for their PHI Access	08
4.2	Storage and Removal of PII	08
05	Automatic Log-Off	08
06	Audit Controls	08
07	Authorization	09
7.1	Access Prevention	10
7.2	Privilege and Access	10
08	Cyber Security Product Upgrades	10
09	Health Data De-Identification	11
10	Data Backup and Disaster Recovery	11
11	Health Data Integrity and Authenticity	11

Table of Contents

12	Malware Detection/Protection	11
	12.1 Other Compensation/Protection Controls	11
	12.2 Firewall Implementation	12
13	Connectivity Capabilities	12
14	Personal Authentication	12
15	Roadmap for Third Party Components in Device Life Cycle ...	13
16	System and Application Hardening	13
17	Health Data Storage Confidentiality	14
18	Transmission Confidentiality	14
19	Security Program Integration	15
	19.1 Risk Management	15
20	Secure Decommissioning	15

01 Purpose

This Security Operations Manual (SOM) provides information that Stryker's customers need to know in order to integrate a specific Stryker device or health IT solution into a customer's IT network environment in a secured manner.

It also supports customer's ability to perform risk management, to identify configurable security controls, and to better protect their systems.

02 Definitions

API – Application Programming Interface

An interface for computing that defines interactions between multiple software intermediaries.

COTS – Commercial off-the-shelf

Software (or any other item) that is sold as a packaged solution which is then adapted to satisfy the needs of the organization purchasing the COTS. Some medical devices utilize COTS software in addition to or instead of software developed by the manufacturer.

Refer- third-party software.

Customer

The individual or organization responsible for procurement and operation of the device. See Owner and Operator.

Device

The item being integrated or used for a healthcare purpose. A Medical Device or other health IT product may be referred to as a Device or a Product in this document.

DICOM (Digital Imaging and Communications in Medicine)

Standard developed by NEMA and the American College of Radiology, used worldwide to store, exchange, and transmit medical images.

FDA – U.S. Food and Drug Administration

A federal agency of the United States' Department of Health and Human Services.

Refer www.fda.gov

HDO – Healthcare Delivery Organization

Also "Health Delivery Organization," an organization or group of organizations that are involved with the delivery of healthcare services. A hospital is an HDO. If an HDO purchases and operates a Stryker device, the HDO is also the Customer, Owner, and Operator per the definitions of those terms.

IEC – International Electrotechnical Commission

A global organization whose work underpins quality infrastructure and international trade in electronic goods. IEC publishes thousands of international standards, including documents related to medical device software (for example, IEC 62304).

Refer www.iec.ch.

IFU – Instructions for Use

Information provided by the manufacturer in document or electronic form, informing the user of a device's intended purpose and proper use and of any precautions to be taken.

ISO – International Organization for Standardization

An international standard-setting body that promotes proprietary, industrial, and commercial standards, and publishes standards relevant for information technology, privacy, and security (for example, ISO/IEC 27034).

Refer www.iso.org

Manufacturer

The entity (Stryker) that builds the device and sells it to the customer.

MDR – European Union (EU) Medical Device Regulation of 2017

The European Union regulation concerning medical devices.

Refer https://ec.europa.eu/health/md_sector/overview_en

Medical Device

See the following sources if a precise definition is required: FDA, MDR (EU) 2017/745, ISO 14971:2007.

NEMA – National Electrical Manufacturers Association

Refer www.nema.org

NIST - National Institute of Standards and Technology

A physical sciences laboratory and non-regulatory agency of the United States Department of Commerce. NIST has published comprehensive standards for the selection, implementation, and risk management of security and privacy controls (e.g., NIST SP 800-53).

Refer www.nist.gov.

Operator

The person(s) using the device for its intended purpose. This term may also sometimes refer to the person or organization responsible for procuring the device (owner, customer).

OSS – Open-Source Software

Third party software licensed under an OSS license, in which the copyright holder grants users the rights to use, study, change, and distribute the software to anyone and for any purpose as long as the license terms are adhered to.

Owner

See Operator and Customer.

PHI - Protected Health Information

Individually identifiable health information (IIHI) that is transmitted by electronic media; maintained in electronic media; or transmitted, or maintained, in any other form or medium (source: extracted from 45 CFR Section 160). Note: This is a subset of PII.

PII - Personally Identifiable Information

Any information about an individual maintained by an agency, including the following:

- Any information that can be used to distinguish or trace an individual's identity.
- Any other information that is linked or linkable to an individual, such as medical, educational, financial, and employment information (source: from NIST SP 800-122).

Product

See Device.

SaMD - Software as a Medical Device

Software intended to be used for one or more medical purposes that perform these purposes without being part of a hardware medical device (source: from International Medical Device Regulators Forum).

SBoM – Software Bill of Materials

For a specific device, a listing of all software components that are incorporated into the final product. The SBOM may be used to assist with operational security planning by the HDO.

SOM - Security Operations Manual

A product-specific guide to the secure integration of a product into a customer IT network (this document).

Third-party software

Third party software is software not developed by Stryker, and for which Stryker otherwise does not have complete ownership. See COTS and OSS.

User

See Operator.

03 Product Description

This Security Operations Manual (SOM) provides information that Stryker's customers need to know in order to integrate a specific Stryker device or health IT solution into a customer's IT network environment in a secured manner.

It also supports customer's ability to perform risk management, to identify configurable security controls, and to better protect their systems.

Manufacturer Name	stryker®
Stryker Division	Stryker Global Technology Center
Address	Stryker Global Technology Center Private Limited 10th Floor, Vatika Business Park, Block Two, Sector-49 ,Sohna Road, Gurgaon 122002, Haryana, India
Device Description	Ortho Slice 3D Knee Planning software is used to create a pre-operative planning for a knee replacement surgery where Stryker's Triathlon knee implant is used. The Ortho Slice 3D Knee Planning software is intended to provide a surgeon facing, easy to use knee planning software, that uses the patient's CT scans to visualize the disease condition of Knee in three-dimension and enable effective decision making for the surgeons before they even go into the operating room on the day of the surgery.
Device Model, Version	6007-670-000 V1.0 (Further digits for minor fixes controlled internally)
Manufacturer Contact Information	Manufacturer: Stryker Global Technology Center Private Limited 10th Floor, Vatika Business Park, Block Two, Sector-49, Sohna Road, Gurgaon 122002, Haryana, India Distributed By: Stryker Japan K.K. 2-6-1, Koraku, Bunkyo-ku,Tokyo, 112-004, Japan t/f: 03-6894-0000 Additional information and contact links are available on Stryker's Product Security webpage, https://www.stryker.com/us/en/about/governance/cyber-security.html .

Table 1.1 Product Description

3.1 Device and Manufacturer Identification

Device:

Ortho Slice 3D Knee Planning Software

Manufacturer:

Stryker Global Technology Center Private Limited

10th Floor, Vatika Business Park
Block Two, Sector-49, Sohna Road,
Gurgaon 122002, Haryana, India

3.2 Device Intended Use:

The Ortho Slice 3D Knee Planning software is intended to provide a surgeon facing, easy to use knee planning software, that uses the patient's CT scans to visualize the disease condition of the knee in three dimension and enable effective decision making for the surgeons before they even go into the operating room on the day of the surgery.

Functionality Includes:

- Auto-segmentation and landmark identification (manual modification possible)
- Effortless implant planning
- Saves the planned report for quick reference, the surgeon has the control and can choose how to interpret and use the results from the pre-operative planning.

Contraindication:

The surgeon needs to determine whether the patient's conditions are appropriate for this kind of procedure. The patients who have other type of metallic implants at or near Knee, Ankle, or Hip joints, which can create artifacts/noise on the CT Scan, are against the use of this system. In addition, some patients with advanced osteoporosis or deformities would be contraindicated like fused bone in the femur and tibia at the knee region. Reads only CT data.

Refer IFU and User Manual for more details.

3.3 Vulnerability Intake and Monitoring:

When Stryker obtains vulnerability information through surveillance or other sources, an assessment of the vulnerability's exploitability and impact is conducted. Based upon this assessment Stryker determines if further actions are required like providing security updates and/or providing communication to the customer in a timely manner. Vulnerability information may also be requested from Stryker at any time.


Any potential security vulnerabilities customer may become aware of due to Ortho Slice 3D Knee Planning software must be communicated to Stryker customer care and the same will be handled through the post market complaints management process to do the assessment and required actions including any updates needed for the customers.

3.4 System Characterization and System Assets

Ortho Slice 3D Knee Planning software allow surgeon to create preoperative planning before proceeding to the surgery (Conventional or Navigated). This application allows to import/load Patients DICOM CT images from external storage devices such as USB, Hard disk and network drive. This application will not allow user to transfer the patient data to any other external or connected system to process further. All the patient data is encrypted and stored locally under the logged user folder.

3.5 System Security Context and Intended Environment

Figure 1: System Security

While there is no specific requirement for Ortho Slice 3D Knee Planning software to be fully functional other than a usual windows environment, however Stryker recommends the user to follow some of the best practice security standards in order to run the application in a safe and secure environment as follows:

Devices operating in the intended use environment should consider that their IT infrastructure must follow different risk management approaches associated with their networks. HDO or customer must adopt a risk management process adhering to general cybersecurity best practices to maintain the healthcare provider's overall security status and their secure environment, as follows:

- Good physical security to prevent unauthorized physical access to Ortho Slice 3D Knee Planning software application.
- Access control measures (for example, role based) to ensure only authenticated and authorized personnel are allowed access to network elements, stored information, services and applications.
- General patch management practices that ensure timely security patch updates.
- Malware protection to prevent unauthorized code execution.
- Security awareness training.

3.6 Setup of the SaMD (Software as a Medical Device) Ortho Slice 3D Knee Planning

The Ortho Slice 3D Knee Planning software operates on Window-based PC

Specifications	Description
Operating System	Windows 10 version
CPU	Intel Core i7-4770 or higher
RAM	32 GB or more
Disk space	12 GB for program installation, 500GB or more for patient data storage
Monitor resolution	1920*1080 pixels or above
Graphics	Dedicated 4GB or higher memory (for example, Nvidia Quadro RTX3000 or NVIDIA T600)
Communication port	USB 2.0 or higher
Optional mouse with left and right button, and scroll	

Note:

- For cases where the patient data is loaded from local network drive, ensure network connectivity before starting the system. If there is a disruption in network connectivity, loading the data may take longer time than usual. If there is a disruption in network while loading the case, the software will notify the user.
- Information for compatibility with other devices. External hard drive to transfer CT data should be compatible with USB 2.0 or higher.
- For other Security related requirements please refer this manual

04 Management of PII and PHI

Ortho Slice 3D Knee Planning software does not process PII /PHI outside the surgeon's laptop. The HDO has full control of it is laptop including the Ortho Slice 3D Knee Planning software and responsibility for any PII /PHI there. The software will only display the PII /PHI information from the DICOM CT data or entered by the Surgeon and will not process it outside of the surgeon's laptop boundary. The HDO and User has the full control of the PHI and PII data in the Laptop and if any data to be removed based on the HDO data retention Policies, the HDO/ User can take support from Stryker to remove the data not needed to be retained.

4.1 Handling of Patient Requests for their PHI Access

“Refer to Section Management of PII and PHI” above. User/ HDO has full control on patient data kept on laptop. Application do not have any additional functionality to provide access for patient requests. The HDO need to take the Report Output and share with Patients in cases where the Patient request his Personal Health Information as per the HDO process

4.2 Storage and Removal of PII

Stryker does not process outside the surgeon's laptop. The HDO has full control of it is laptop including the Ortho Slice 3D Knee Planning software and responsibility for any PII there.

PII is embedded in the input DICOM files provided by the surgeon to the software. The software uses the PII data to display on its GUI and final planning PDF contains the PII information from DICOM files.

PII data will be maintained in volatile memory within the user laptop and can be exported to a PDF file. PII data is not transferred to any other system. User has the provision to anonymize the Patient Name on the Software GUI for the purpose of presenting the data for any external stake holders for training or other needs to support any privacy requirements for the HDO.

For backup and restore of PII data, please reach out to Stryker customer care for support if needed.

05 Automatic Log-Off

Customers are advised to configure windows OS to automatically lock the screen after a period of idle time as per the HDO IT policies.

Application also has the ability to lock the screen after inactivity for configurable timeout. User can configure the inactivity timeout. For details, please refer User Manual Section To configure System Settings.

06 Audit Controls

Stryker uses strong protection mechanism to protect the audit logs from getting tampered by any unauthorized party and hence does not require any extra steps from the users. Audit logs are encrypted using 256-bit AES encryption to avoid any tampering of the information. Decryption of the audit log is handled by Stryker on request from authorities. User cannot edit or alter the audit logs.

Ortho Slice 3D Knee Planning Software captures the following type of audit events:

- Creation/modification events of patient PII data (No PII data stored)
- Import of DICOM data from removable media
- Application Programming Interface (API) and similar activity – Used for Printing and PDF view.
- Marked data with time stamp information to enable it to be selected for deletion based on when it was acquired or stored. Work step visited, save / update operation on patient data. Application data - workflow and features executed.

Audit Logs Format is: <timestamp>,<user>,<component>,<Feature/Module>,<Action>

It is possible to export the logs via physical media considering the physical media like USB etc. to be secure. But it is recommended for the users to keep their physical media secure and updated against the latest threats.

Below are some of the safety measures that can be implemented to secure physical media like USB drive.

Do not plug a USB drive into an unknown computer

Do not plug the USB into any computer without verifying the identity and safety of that computer system as the system may pose a potential security threat to your physical device.

Take advantage of security features

Use passwords and encryption on your USB drive to protect your data, and make sure that you have the information backed up in case your drive is lost.

Disable Autorun

The Autorun feature causes removable media such as CDs, DVDs, and USB drives to open automatically when they are inserted into a drive. Disabling Autorun prevents malicious code on an infected USB drive from opening automatically.

Use and maintain security software and keep all software up to date

Use a firewall, antivirus software, and anti-spyware software to make your computer less vulnerable to attacks and keep the virus definitions current. Also, keep the software on your computer up to date by applying any necessary patches.

07 Authorization

Ortho Slice 3D Knee Planning software requires a valid license in order to be fully functional and running., which needs to be obtained from Stryker only. Apart from requiring a valid license, the application user must leverage windows authentication and authorization mechanism for user access to software and data. So, Stryker best recommends the user to setup proper authorization of users on their laptops as discussed below.

Authorization in system security is the process of giving the user permission to access a specific resource or function. In secure environments, authorization must always follow authentication. Users should first prove that their identities are genuine before an organization's administrators grant them access to the requested resources. So proper authorization must be implemented at system level to harden the security. Different approaches to authorization may include

Role-Based Access Control (RBAC)

Users are identified as being in a role that stipulates what privileges they have. Additionally, their user ID would restrict what data they have access to.

Access Control Lists (ACL)

An ACL specifies which users have access to particular resources. For instance, if a user wants to access a specific file or folder, their username or details should be mentioned in the ACL in order to be able to access certain data.

7.1 Access Prevention

Ortho Slice 3D Knee Planning software does not have any built-in access prevention features enabled in the Ortho Slice 3D Knee Planning software and leverages windows access prevention mechanism. It is recommended to customers to have proper access control measures as discussed in the below sections.

Taking steps to prevent unauthorized access to the system and its software components is important for a wide number of reasons, including preventing others from installing spyware and deleting your important files, or even creating viruses. By making changes to your computer to prevent unauthorized access, you are also protecting your personal privacy. Here are some steps to take to properly secure your computer and prevent others from accessing or modifying your application data:

- Set up password protection for user authentication: Password protection must be enabled at system level so that any unauthorized user cannot access the system. Password should be set in such a way that it must not be easy to guess. Also, strong password policy must be implemented to enhance the overall security of the system
- Install antivirus software or a spyware protection program: A good antivirus or spyware protection program must be installed on the system. These programs are used to detect any malicious actions or programs that might be used as a threat for the system or installed application. Lastly these antivirus programs must be regularly updated so that it can protect the system and the installed applications from latest security threats.
- Restrict the access to your system only to a limited number of trusted peoples. This can help the installed application to be accessed only by an authorized individual.

7.2 Privilege and Access

Stryker recommends the laptop administrator to create separate users with appropriate privileges for access to Ortho Slice 3D Knee Planning software on the same laptop. Privilege and access to the Ortho Slice 3D Knee Planning software must be restricted such that any user of the application can only use it within its intended use and any other functionality outside of the scope of the application is restricted as much as it can be. Users can maintain their own data on same laptop without access to other user's data by setting different log in access.

08 Cyber Security Product Upgrades

The application does not have any updates installation policy implemented. Hence the users will not get any online updates. If any potential vulnerabilities are identified by Stryker which require an update at the customer site, a new version of the software will be released, and customers will be informed about the action to be taken at their end.

It is HDO's responsibility to update the latest patches for their operating system, their third-party components (if any) and other applications like Virus Protection software's/anti-malware software's, firewalls etc. Timely to ensure the security and protection of the system.

Ortho Slice 3D Knee Planning software does not contain any malware protection embedded in it. Hence, users are advice to install and anti-malware software on laptop.

09 Health Data De-Identification

Ortho Slice 3D Knee Planning software anonymizes the data at runtime but does not delete or remove them. This can be done via GUI interface of the application itself.

Refer User Manual Section Additional tools for more details.

10 Data Backup and Disaster Recovery

The application does not contain any online or offline mode of data backup or its recovery. So, the users are expected to have their own copy of data backup possibly in any physical media or via some online storage methods.

11 Health Data Integrity and Authenticity

No user actions are needed since any health data or other sensitive data stored on the system is encrypted using strong 256-bit AES encryption algorithm by the application itself to preserve the data integrity. The application properly checks the integrity of the data before loading them.

12 Malware Detection/Protection

The standalone Ortho Slice 3D Knee Planning software by default does not contains any malware detection functionality and requires the user to have some malware detection in place. As, the malware detection is crucial with malware's prevalence because it functions as an early warning system for the computer secure regarding malware and cyberattacks. It keeps hackers out of the computer and prevents the information from getting compromised. This involves the process of scanning the computer and files to detect malware.

To protect against the malwares below points are recommended:

- Install a good malware detection program on the system
- Keep your computer and software updated
- Use a non-administrator account whenever possible
- Think twice before clicking links or downloading anything
- Be careful about opening email attachments or images
- Do not trust pop-up windows that ask you to download software
- Limit your file sharing

12.1 Other Compensation/Protection Controls

The Ortho Slice 3D Knee Planning Software application contains several protection mechanisms by its design like the application requires a valid license in order to work properly, the logs are encrypted using private keys and the data de-identification is also done which anonymizes the data at runtime. Besides these projections in place, Stryker recommend the users to apply certain other safety measure as below:

- The application should be allowed to be run only as authorized individual enabled using built in mechanism of windows OS.
- Proper Application whitelisting should be done on all security agents running on the device so that the

Ortho Slice 3D Knee Planning software is not flagged as malicious in any case.

- Any third-party components installed in the system must be properly updated.
- Regular Antivirus scan should be done in order to eliminate any possible threats.
- Application Logs shall be audited for any errors or proper functioning of the application.
- Regular windows updates and patches should be installed.

12.2 Firewall Implementation

In order to safely use the application its important the below listed firewall rules be strictly followed and respected while using the application:

Ensure that you have enabled the firewall. If not, please follow the below steps on your Windows computer

1. From the Start menu, click Control Panel
2. Select System and Security. Click Windows Defender Firewall
3. Click Turn Windows Firewall on or off
4. Select Turn on Windows Firewall for domain home/work (private) or public network

Other things to consider:

- It is imperative that logs/warnings get the attention of the user. Any anomaly must then be resolved after its addressed.
- Ensure that ports necessary for function are accessible only to authorized clients of the application.
- In case a malware was reported on the system, ensure a proper sweep scan has been initiated and the removal of the malware was successful before resuming normal operations.

Firewall helps in preventing network access to devices. If properly used and configured it can lead to protected and reliable accessibility. It can help in prevention of unauthorized access and network connections against external threats, IP spoofing & routing attacks and malicious packets.

13 Connectivity Capabilities

The Ortho Slice 3D Knee Planning software application by default does not require any network connectivity or even wireless connectivity for its operation. The input CT data can be loaded from a network drive and the HDO IT to ensure adequate security controls to protect the network drive connected to the Laptop where Ortho Slice 3D Knee Planning software is installed. Also, it does not apply any restrictions in place when the physical media needs to be inserted in the system for data backup or data transmission.

So, it is advised for the users to use only secure and updated physical media in case it is required. To secure the USB or any physical media its best to properly scan them using a good antivirus program.

14 Personal Authentication

Ortho Slice 3D Knee Planning software does not provide any authentication mechanism apart from requiring a valid license which is unique for the system in order to be operational.

Stryker recommend the user to implement the secure windows authentication using strong password-based authentication on their laptop. These passwords should be strong enough which is not easy to guess. Also, it must contain alphanumeric characters along with special characters to ensure best security practice.

For proper user management, the system should be configured for the below points.

- The authentication system should be done via password-based login or integrated windows-based authentication
- After a few unsuccessful attempts account must lockout
- Passwords must be changed after a regular interval of time
- Default password or easy guess password must not be accepted by the system. In other words, the password should meet the password complexity policy
- The system must be configured in such a way as to lock if it is left idle after a reasonable period
- Physical security must also be implemented to manage access to the system

15 Roadmap for Third Party Components in Device Life Cycle

Stryker has evaluated third -party components as per the requirement identified in IEC 62304 and adequate actions are implemented in application.

Stryker will be evaluating high-risk third-party components periodically and communicate to customers for any updates required during the product lifecycle.

16 System and Application Hardening

Stryker had performed the application security testing and security code review of Ortho Slice 3D Knee Planning software. Ortho Slice 3D Knee Planning software is hardened by eliminating any vulnerability or flaw, which can lead to security issue.

Systems hardening is a collection of tools, techniques, and best practices to reduce vulnerability in the application, systems, and other areas. The goal of systems hardening is to reduce security risk by eliminating potential attack vectors and condensing the system's attack surface. By removing superfluous programs, accounts functions, applications, ports, permissions, access, so on. attackers and malware have fewer opportunities to gain a foothold within the IT ecosystem. Systems hardening demands a methodical approach to audit, identify, close, and control potential security vulnerabilities. The type of hardening carried out depends on the risks in the existing technology, the resources that are available, and the priority for making fixes.

Stryker recommends to customers to keep below key points while implementing the system hardening.

Audit your existing systems

Carry out a comprehensive audit of the existing technology. Use penetration testing, vulnerability scanning, configuration management, and other security auditing tools to find flaws in the system where the application is installed and prioritize fixes

Create a strategy for systems hardening

There is no need to harden all the systems at once. Instead create a strategy and plan based on risks identified within the technology ecosystem and use a phased approach to remediate the biggest flaws.

Patch vulnerabilities immediately

Ensure to have an automated and comprehensive vulnerability identification and patching system in place.

Network hardening

Ensure the firewall is properly configured and that all rules are regularly audited; secure remote access points and users; block any unused or unneeded open network ports; disable and remove unnecessary protocols and services; implement access lists; encrypt network traffic.

Also refer to Section 13

Operating system hardening

Apply OS updates, service packs, and patches automatically; remove unnecessary drivers, file sharing, libraries, software, services, and functionality; encrypt local storage; tighten registry and other systems permissions; log all activity, errors, and warnings; implement privileged user controls.

Eliminate unnecessary accounts and privileges

Enforce least privilege by removing unnecessary accounts (such as orphaned accounts and unused accounts) and privileges throughout the IT infrastructure.

Anti-Malware installation

The system running Ortho Slice 3D Knee Planning software should have proper anti-malware software installed with latest updates.

17 Health Data Storage Confidentiality

The data at rest is encrypted using a strong encryption mechanism implemented within the application which safeguards the sensitive medical data from prying eyes.

18 Transmission Confidentiality

Ortho Slice 3D Knee Planning software does not transmit the data over network or internet. End user inputs data into the software using removable media or from local system. The pdf output can be stored or exported on the removable media.

Below are some of the guidelines to the user to be followed while managing data confidentiality:

Manage data access

Controlling confidentiality is, in large part, about controlling who has access to data. Ensuring that access is only authorized and granted to those who have a “need to know” goes a long way in limiting unnecessary exposure. Users must also authenticate their access with strong passwords and, where practical, two-factor authentication. Periodically review access lists and promptly revoke access when it is no longer necessary.

Physically secure devices

Controlling access to data includes controlling access of all kinds, both digital and physical. Protect devices from misuse or theft by storing them in locked areas. Never leave devices or sensitive documents unattended in public locations.

Securely dispose of data

When data is no longer necessary for any-related purposes, it must be disposed of appropriately.

Manage data acquisition

When collecting sensitive data, be conscious of how much data is actually needed and carefully consider privacy and confidentiality in the acquisition process. Avoid acquiring sensitive data unless necessary; one of the best ways to reduce confidentiality risk is to reduce the amount of sensitive data being collected in the first place.

Manage data utilization

Confidentiality risk can be further reduced by using sensitive data only as approved and as necessary. Misusing sensitive data violates the privacy and confidentiality of that data and of the individuals or groups the data represents.

Manage laptop

Computer management is a broad topic that includes many essential security practices. By protecting devices, you can also protect the data they contain. Follow basic cybersecurity hygiene by using anti-virus software, routinely patching software, whitelisting applications, using device passcodes, suspending inactive sessions, enabling firewalls, and using whole disk encryption.

19 Security Program Integration

Ortho Slice 3D Knee Planning software is a standalone software installed on user laptop. The secure practices for the software are covered under Section 8 of this document.

Stryker has implemented incident response programs as part of complaint handling process for the product. The product is intended for Japan market only and user can call t/f: 03-6894-0000 to inform Stryker if any potential threat is identified in the software. Then, Stryker team will take the right steps to help the customer to deal with the situation. Stryker has done extensive security testing of the Ortho Slice 3D Knee Planning software application and implemented adequate actions to ensure protection from external threats. However, beyond this security measures in place, it is advised for the users to take a step ahead and follow some of the below guidelines to ensure better security postures.

- Do not plug any unknown physical media like USB etc. If it is required to plug in to the device, it must be scanned thoroughly using a strong anti-malware program.
- System must be scanned on a regular basis for any potential threats using Anti-malware and/or anti-virus softwares

19.1 Risk Management

Stryker integrates cyber security risk management into its overall program for health and safety risk management. Both security and safety risk assessments were conducted for this device per guidelines in compliance with EN/ISO 14971 and Stryker Product Security procedures. Additionally, Stryker has a robust post-market security risk management process that monitors the ongoing security posture of this device and addresses any security incidents that might arise.

20 Secure Decommissioning

Please reach out to Stryker Customer Care for secured decommissioning.

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Stryker Corporation or its divisions or other corporate affiliated entities own, use or have applied for the following trademarks or service marks: Ortho Slice 3D Knee Planning, Triathlon. All other trademarks are trademarks of their respective owners or holders.